



Growing X-ray Quality Crystals



*Paul D. Boyle Department of Chemistry
University of Western Ontario*

and

*Louise N. Dawe
Chemistry and Biochemistry Department
Wilfrid Laurier University*

For a more in-depth discussion:

- [Paul Boyle's CCCW22 slides](#)
- <http://xray.chem.uwo.ca/Guides.html>
- <http://www.nottingham.ac.uk/~pczajb2/growcrys.htm>
- <http://www.cryst.chem.uu.nl/lutz/growing/growing.html>
- <http://xrayweb.chem.ou.edu/notes/xtalgrow.html>
- [Scientists grow diamonds from scratch in 15 minutes](#)

Important Caveats!

- X-ray crystallography **does not** determine your compound's chemical composition
- It is a technique for determining atomic connectivity of a single crystal.
- Is your single crystal representative of the bulk composition?
- Samples submitted for X-ray crystallographic analysis need to have other characterization data (e.g. NMR, IR, MS, ESR, elemental analysis, etc.) to establish the chemical formula
- Do not make assumptions or jumps in logic that your crystal structure is actually representative of your major reaction product.

What Are Good Crystals and Why We Need Them




A “good” crystal:

- is 0.1 - 0.3 mm in at least two of its dimensions
- exhibits a high degree of internal order as evidenced by the presence of an X-ray diffraction pattern
- very often, but not always, shows regular faces and edges




We need “good” crystals because:

- Quality of sample is characterized by maximum diffraction angle (θ); also expressed in “resolution” (Å)
- The larger the max. diffraction angle, the higher the resolution and the greater number of data (which are necessary to adequately model the structure)
- Discerning individual atomic positions requires data resolution which is higher than chemically significant distances (e.g C=O \sim 1.2 Å)

The Effect of Limiting the Resolution of the Data

<p>The final structure</p>	<p>Electron density map using all available data ($\theta_{\max} = 32.35^\circ$)</p> <p>Resolution = 0.66 Å</p> <p>All atomic positions are easily resolved</p>	<p>Limited $\theta_{\max} = 25.0^\circ$</p> <p>Resolution = 0.84 Å</p> <p>Peaks are beginning to flatten out</p> <p>Atomic positions are still easily resolvable</p> <p>IUCr recommended minimum resolution</p>
		

The Effect of Limiting the Resolution of the Data

<p>Limited $\theta_{\max} = 19.47^\circ$</p> <p>Resolution = 1.5 Å</p> <p>Peaks start to “melt” into each other</p> <p>Individual atomic positions are still resolvable</p>	<p>Limited $\theta_{\max} = 14.48^\circ$</p> <p>Resolution = 2.0 Å</p> <p>Metal position resolvable</p> <p>Only gross shape of organic ligand evident</p> <p>Peak positions for ligand have shifted away from true atomic positions</p>	<p>Limited $\theta_{\max} = 7.18^\circ$</p> <p>Resolution = 4.0 Å</p> <p>Metal position is only barely above background</p> <p>No trace of ligand</p>
		

The Right Attitude Toward Crystal Growing for X-ray Analysis

- Growing X-ray quality crystals requires care and attention to detail
- Treat it like its own miniature research project
- If you have spent weeks or months doing your synthesis, why assume finding the correct crystallization conditions will just take a few hours?
- Do not try to skimp on the amount of material when growing crystals
- Purify your compound (using conventional crystallization and/or other purification steps)
- Consider the empirically established physical properties of your compound – sensitivities, thermal stability, etc.
- Develop a solubility profile of your compound
- Use CLEAN glassware as crystal growing vessels
- Set up crystal growing attempts in parallel utilizing different conditions

Special Considerations for Glassware

- Before setting up a crystal growing attempt think about how the crystals will be handled
- Crystals will need to be extracted from the vessel without damage
- Therefore, pick a suitable crystal growing vessel

BAD: Round bottom flasks of any size

BAD: Small aperture vials (too small for spatula)

BAD: Screw top vials (shoulder causes difficulty)

GOOD: NMR tubes

GOOD: Small test tubes

GOOD: Vials without shoulders

Factors Affecting Crystallization

- Solvent – moderate solubility is best. Supersaturation leads to sudden precipitation and smaller crystal size
- Nucleation – fewer nucleation sites are better. Too many nucleation sites (i.e. dust, hairs, etc.) lower the average crystal size
- Mechanics – mechanical disturbances are bad.
- Time – faster crystallization is not as good as slow crystallization. Faster crystallization results in a higher chance of lower quality crystals

Crystal Growing Techniques

Slow Evaporation: simplest to set up.



- Has drawbacks: solute can “oil out”, crystals stick to sides of vessel making them difficult to extract from vessel without breaking them.


Slow Cooling: Soluble when hot, insoluble when cool.



- Use Dewar to slow the cooling process.



Variations: use binary or tertiary solvent mixtures.

- Use solvent with similar boiling points and other properties.
- Document the percentages of each solvent component!

Vapour Diffusion	Solvent Diffusion	Reactant Diffusion
<p>dissolve solute in solvent (green)</p> <p>precipitating “anti- solvent” (pink)</p> <p>pink should be more volatile than green</p> <p>Do not let sides of small vessel touch vertical surface of outer vessel (prevent capillary action)</p>	<p>Good for milligram amounts</p> <p>Use NMR tube for best results</p> <p>Fill soluble more dense solvent on bottom with your solute.</p> <p>Fill the rest of tube with less dense precipitant solvent E.g. $\text{CH}_2\text{Cl}_2/\text{Et}_2\text{O}$</p>	<p>Set-up similar to solvent diffusion except that reactants are in different layers</p> <p>Good for milligram amounts</p> <p>Good for completely insoluble products which never go back into solution after being formed</p> <p>Consider using a 3rd “middle layer” solvent to mediate the reactant concentrations</p>
 <p>The diagram shows a cross-section of a vapour diffusion apparatus. It consists of an 'Outer vessel (closed)' which is partially filled with a pink 'Anti-solvent'. Inside this is an 'Inner vessel' which is partially filled with a green 'Solution'. The inner vessel is inverted so that its open end is submerged in the anti-solvent. At the top of the inner vessel, there are 'Gaps' between the vessel walls and the stopper. The green solution is shown as a column inside the inner vessel, with some green droplets indicated near the interface with the anti-solvent.</p>	 <p>The diagram shows a cross-section of a solvent diffusion apparatus, which is a vertical tube. The bottom portion of the tube is filled with an orange liquid, representing the solvent containing the solute. The top portion of the tube is filled with a blue liquid, representing the precipitant solvent. A horizontal interface is visible between the two liquids, with small purple droplets shown at the interface, indicating the formation of a product.</p>	

Slow Cooling	Convection (Principles)	Convection (Easy Way)
<p>Heat water to just below the boiling point of the solvent you plan to use</p> <p>Create a saturated solution of your compound in boiling solvent</p> <p>Transfer to solution to a test tube and cap</p> <p>Fill Dewar with hot water</p> <p>Place test tube into the Dewar and cap Dewar</p> <p>Let stand several days</p>	<p>Good for insoluble or sparingly soluble compounds</p> <p>Create a thermal gradient in the crystal growing vessel</p> <p>Solvent becomes saturated in “warm” region and deposits material in “cool” region where nucleation and crystal growth can occur.</p>	<p>Local cooling – simple to set up</p> <p>Take flat bottomed crystal growing dish and set up like slow evaporation</p> <p>Place vessel so that one side is against a heat sink, e.g. an outside window (in Winter at least)</p>
<div data-bbox="254 856 662 1299"><p>Capped vial</p><p>Crystals</p><p>Saturated solution</p><p>Thermal reservoir</p></div> <p data-bbox="254 1328 662 1370">Crystal growth by controlled cooling</p>	<p>Cyclic current allows continual replenishing of the solute</p> <p>Velocity of convection current is proportional to the magnitude of thermal gradient.</p> <p>Take care that the gradient isn't too large – too high velocity inhibits crystal growth</p>	<p>Placing crystal growing dish on a cool surface will not cause convection</p>

Hope's Convection Device	Thiele Tube Convection Device
<p>Specialty glassware can be fabricated</p> <p>This device designed by Håkon Hope (J. Appl. Cryst. 1971, 4, 333)</p> <p>Nichrome wire is used for the main [A] and subsidiary [B] heating elements.</p> <p>The solvent is cooled at point [C] using Cu wire or thin plastic tubing as a heat sink.</p>	<p>If you don't want to fabricate a special piece of glassware – improvise!</p> <p>Fill Thiele tube with solvent. Wrap nichrome wire around the bottom side arm and attach to Variac</p> <p>Place solute in small container just below top side entrance</p> <p>Apply heat</p>
 <p>A schematic diagram of Hope's Convection Device. It shows a vertical glass tube with a side arm at the bottom. The main vertical tube has a heating element (nichrome wire) labeled 'A' at the bottom. The side arm has a heating element labeled 'B' and a cooling point labeled 'C' at the top. The diagram is labeled with a small '3' in the top right corner.</p>	 <p>A photograph of a Thiele Tube Convection Device. It is a glass Thiele tube with a side arm at the bottom. The side arm is wrapped with nichrome wire, which is attached to a Variac. The tube is filled with solvent, and a small container with solute is placed just below the top side entrance. The device is shown against a purple background.</p>

Soxhlet Extraction	Sublimation	More comments on Sublimation
<p>Soxhlet extraction is normally used for separations.</p> <p>However, it can be used for crystal growing of thermally stable, sparingly soluble materials</p> <p>Place your solute in the sample thimble</p> <p>Start the solvent refluxing</p> <p>Crystals grow in the solvent reservoir</p>	<p>Good for volatile air- sensitive materials</p> <p>Specialty glassware is available</p> <p>Use minimal heat to sublime slowly</p> <p>Use small amounts of material</p>	<p>Gas to solid phase crystal growth</p> <p>Compound needs to be thermally stable</p> <p>Can be easy to set up – vacuum sealed tube of material placed in oven for several days/weeks</p> <p>Or more complicated – material packed in tube followed by glass wool.</p>
		<p>Place under active or static vacuum and set-up thermal gradient by heating the loaded end of the tube.</p> <p>Place Cu pipe around tube to create thermal gradient.</p>

Chemical Modification	Ionization of Neutral Compounds	Co-Crystallants
<p>For ionic compounds, change the counterion to change the solubility and other characteristics of your compound</p> <p>Ions of similar sizes tend to pack together better</p> <p>Use counterions with rigid geometries e.g. triflate, BPh_4^-, Me_4N^+, $(\text{Ph}_3\text{P})_2\text{N}^+$</p> <p>Tend to disorder: Et_4N^+, Bu_4N^+, BF_4^-, PF_6^-</p> <p>Make sure counterion does not react with your compound!</p>	<p>If your compound is neutral and has proton acceptor or donor groups, consider ionizing the compound</p> <p>The ionic form may take advantage of hydrogen bonding to give better crystals</p> <p>Counterions can be changed to optimize crystal growth</p> <p>This will change your compound, but if you are only interested in confirming a structure, and not in detailed electronic properties, this shouldn't be a problem.</p>	<p>Sometimes two (or more) different compounds “co-crystallize”.</p> <p>Most commonly, this is a solvent molecule.</p> <p>Triphenylphosphine oxide has been used as a co-crystallant for both inorganic and organic compounds.</p> <p>Tetraaryladamantanes have been used as co-crystallant “chaperones”: Krupp, F.; Frey, W.; Richert, C. <i>Angew. Chem. Int. Ed.</i> 2020, 59, 15875-15879.</p> <p>Rami, F.; Stuerzer, T.; Richert, C.; Adam, M. <i>Acta Cryst.</i> 2021, A77, C883</p>

Final Thoughts and Acknowledgments

- The quality and meaningfulness of your X-ray results is directly dependent on the quality of your sample crystal
- You can get information from a bad crystal structure, but it will be difficult to publish and makes for a weaker manuscript
- Take crystal growing as a serious part of your research project – spend the time and effort to be successful.
- There are many solvents and crystal growing techniques available – use them.
- Thank you: Sandy Blake and Doug Powell for allowing PDB to use images from their crystal growing guides.
- Thank you: Clarence Pfluger, Tony Linden, Sandy Blake, Chuck Barnes, and Andrea Sella for sharing crystal growing methods with PDB.